


Simplicity



Precision



Flexibility



Gebauer SL – The Smart Solution
DSAEK / ALK

Precise, Flexible and Simple.

Gebauer Artificial Anterior Chamber

Precise, Flexible and Simple

The "Click-and-Go" mechanism of the **Gebauer SL DSAEK / ALK Set** assures **safe** and **easy** assembly.

The **flap thickness is selectable**, in steps between 300 µm 450 µm. The **flap diameter can be variably adjusted** up to > 10mm.

The **preassembled** disposable single-use sets guarantee **safe** and **repeatable** results. **Various ring sizes** are available to customize the system for each individual patient.



Gebauer Artificial Anterior Chamber
www.gebauer-sl.com

Gebauer SL DSAEK / ALK Set

Sterilizable Gebauer Artificial Anterior Chamber
Single-use disposable Heads
Reliable "Click-and-Go" mechanism
Automated, motor-driven cuts

DSAEK Instruments

Gebauer offers the following instruments:

Anterior Chamber Cannula

DSAEK-Corneamarker 8,0 mm and 9,00 mm

DSAEK-Spatula (Descemet's stripping)

DSAEK-Hook (Descemet rhexis)

DSAEK-Forceps

DSAEK-Glide

DSAEK-Punch with Single-use Trepphine and Punch Block



Gebauer Instruments
www.gebauer-sl.com





Gebauer Punch

Gebauer Punch with Single-use
Trephine and Punch Block.

Please contact us for more information
or to arrange a product demonstration

www.gebauer-sl.com

All technical specifications subject to change

Gebauer Medizintechnik GmbH

Monbachstr. 7/1
75242 Neuhausen
Germany
6000-003

Phone: +49(0)7234/9421-0
Fax: +49(0)7234/9421-20

info@gebauergmbh.de
www.gebauer-sl.de

